

# Model : GP-500P

## Interferential Low Frequency Stimulator 2Patients



<b>Power voltage</b>	AC 220V 50/60Hz
<b>Power consumption</b>	50W
<b>Maximum output voltage</b>	90Vp-p [500ΩLoaded]
<b>Maximum output current</b>	Auto 50mA[500ΩLoaded] Manual 65mA[500ΩLoaded]
<b>Fixed frequency</b>	2500/4000/5000Hz
<b>Flow frequency</b>	2501~2700/4001~4200/5001~5200Hz
<b>Interferential frequency</b>	1~200Hz
<b>Time</b>	1~30minutes
<b>Suction pressure</b>	150,180,250mmHg
<b>Dimension</b>	418(W)*365(D)*190(H)mm
<b>Weight</b>	8kg



**Accessories:** 1 suction cable, 4 suction cups, 1(white/red) pad cable, 4 pads, 2 (short/long) TENS bands and 4 sponges

## Feature:

- **Portable**
- **Applicable for 2 Patients**
- **Various programs included 9 auto modes**
- **Choose one of the 3 fixed frequency waves (2,500Hz/4,000Hz/5,000Hz)**
- **Fine adjustment to adjust the voltage as 0.2 step precisely**
- **Possible to operate red/white poles separately**
- **Safety function** (zero start method , over current-preventing circuit)
- **Mitigating muscle pain**

لوازم جانبی: یک کابل مکش، دو جفت کاپ مکنده، یک کابل (قرمز/سفید) پد، دو جفت پد، دو باند TENS (کوتاه و بلند) و ۴ عدد اسفنج

## ویژگی‌ها:

- ✓ دستگاه قابل حمل
- ✓ دارای دو خروجی و امکان درمان همزمان و مستقل دو بیمار
- ✓ دارای برنامه‌های متنوع و امکان انتخاب ۹ حالت خودکار
- ✓ امکان انتخاب سه فرکانس ثابت (۲۵۰۰، ۴۰۰۰ و ۵۰۰۰ هرتز)
- ✓ امکان تنظیم دقیق ولتاژ با دقت ۰/۲ ولت
- ✓ امکان تنظیم مجزا از هم قطب‌های قرمز/سفید
- ✓ عملکرد ایمن (روش شروع از صفر، مدار جلوگیری از جریان بیش از حد)

## موارد درمان:

- ✓ کاهش دردهای عضلانی

