

# Model : GP-400D

Low Frequency Stimulator  
2Patients / TENS



<b>Power voltage</b>	AC 220V 60Hz
<b>Power consumption</b>	70W
<b>Maximum output voltage</b>	60Vp-p [500ΩLoaded]
<b>Maximum output current</b>	30mA[500ΩLoaded]
<b>Output frequency</b>	1~1000Hz
<b>Time</b>	1~30minutes
<b>Suction pressure</b>	100, 150, 200mmHg
<b>Dimension</b>	400(W)*350(D)*1100(H)mm
<b>Weight</b>	40 kg



**Accessories:** 2 suction cables, 8 suction cups and 8 sponges

**Sale separately:** 1(white/red) pad cables, 4 pads, 2 (short/long) TENS bands

## Feature:

- Applicable for **2 Patients**
- Various programs included **4 auto modes**
- **16 manual modes** freely selected by the user
- Fine adjustment to adjust **the voltage as 0.2 step** precisely
- Possible to operate **red/white poles** separately
- **Easy to monitor** the status of stimulation **at a glance** through the LCD
- **Easy to move** with 4 wheels embedded under the device
- **Safety function** (zero start method , over current-preventing circuit)
- **Alleviation of muscular pain** by stimulating the percutaneous layer of the skin or scleroderma

لوازم جانبی: دو کابل مکش به همراه چهار جفت کاپ مکند و ۸ عدد اسفنج  
فروش جداگانه: یک کابل (قرمز/سفید) پد به همراه دو جفت پد و دو باند (کوتاه/ بلند) TENS



## ویژگی‌ها:

- ✓ دارای دو خروجی و امکان درمان همزمان و مستقل دو بیمار
- ✓ دارای برنامه‌های متنوع و امکان انتخاب ۴ حالت خودکار
- ✓ امکان انتخاب ۱۶ حالت دستی
- ✓ امکان تنظیم دقیق ولتاژ با دقت ۰/۲ ولت
- ✓ امکان تنظیم مجزا از هم قطب‌های قرمز/سفید
- ✓ پایش آسان و سریع وضعیت برنامه انتخابی از طریق صفحه LCD
- ✓ جابجایی آسان با چهار چرخ تعبیه شده در زیر دستگاه
- ✓ عملکرد ایمن (روش شروع از صفر، مدار جلوگیری از جریان بیش از حد)

## موارد درمان:

- ✓ کاهش دردهای عضلانی با تحریک لایه پوستی یا اسکلرودرمی