

# Model : GP-400P

## Low Frequency Stimulator 2Patients



<b>Power voltage</b>	AC 220V 50/60Hz
<b>Power consumption</b>	40W
<b>Maximum output voltage</b>	30Vp-p [500ΩLoaded]
<b>Maximum output current</b>	8mA[500ΩLoaded]
<b>Output frequency</b>	1~1000Hz
<b>Time</b>	1~30minutes
<b>Suction pressure</b>	120,160,225mmHg
<b>Dimension</b>	418(W)*365(D)*190(H)mm
<b>Weight</b>	8kg



**Accessories:** 1 suction cable, 4 suction cups, 1(white/red) pad cable with 4 pads, 2 (short/long) TENS band and 4 sponges

## Feature:

- **Portable**
- **Applicable for 2 Patients**
- **Various programs included 4 auto modes**
- **16 manual modes** freely selected by the user
- Fine adjustment to adjust **the voltage as 0.2 step** precisely
- Possible to operate **red/white poles** separately
- **Easy to monitor** the status of stimulation **at a glance** through the LCD
- **Safety function** (zero start method , over current-preventing circuit)
- **Alleviation of muscular pain** by stimulating the percutaneous layer of the skin or scleroderma

**لوازم جانبی:** یک کابل مکش، ۴ کاپ مکنده، یک کابل (قرمز/سفید) پد، ۴ پد، دو باند TENS (کوتاه و بلند) و ۴ عدد اسفنج

### ویژگی‌ها:



- ✓ دستگاه قابل حمل
- ✓ دارای دو خروجی و امکان درمان همزمان و مستقل دو بیمار
- ✓ دارای برنامه‌های متنوع و امکان انتخاب ۴ حالت خودکار
- ✓ امکان انتخاب ۱۶ حالت دستی
- ✓ امکان تنظیم دقیق ولتاژ با دقت ۰/۲ ولت
- ✓ امکان تنظیم مجزا از هم قطب‌های قرمز/سفید
- ✓ پایش آسان و سریع وضعیت برنامه انتخابی از طریق صفحه LCD
- ✓ عملکرد ایمن (روش شروع از صفر، مدار جلوگیری از جریان بیش از حد)

### موارد درمان:

- ✓ کاهش دردهای عضلانی با تحریک لایه پوستی یا اسکلرودرمی